

STEP 4. Get result

The blood glucose result will be stored in the memory automatically.



Warning:

1. Please do not change your treatment based on the result without first consulting your healthcare professional.
2. Turn the meter off by removing the test strip. Discard the used test strip and lancet carefully according to your local regulations.


VIEWING THE METER MEMORY

Viewing Results

Your Meter stores the 1000 most recent blood glucose test results with date and time in its memory. You can review the test results with these easy steps.

STEP 1. Enter the memory mode


With the meter turned off, press Main button twice. The first test result will appear, indicating that you are in the memory mode. If you continue to press Main button, you can then review the last 1000 tests in the memory.



When using the meter for the first time "----" will appear, indicating that there are no test results in the memory.

STEP 2. Recall test results

After the last test result, the most recent test result appears with the date and time. Press Main button to review your last 1000 test results in order. When the memory is full, the oldest result will be deleted when the newest is added.



STEP 3. Exit the memory mode

Press the Main button for three seconds to turn off the meter or after the last result, the display will show "End" and it will turn off automatically.




NOTE

The results of control solution tests are **NOT** stored in the memory. Stored results are blood glucose results only.

Reviewing Blood Glucose Day Average Results

1. Press and release Main button.



When **M** appears on the display, keep pressing Main button for 3 seconds until the flashing "R US" appears. Release Main button and then your 7-day average result measured in general mode will appear on the display.

2. Press Main button to review 14-, 21-, 28-, 60- and 90- day average results stored in each measuring mode in the order of Gen, AC, then PC.

3. Exit the meter memory.

Keep pressing the Main button and the meter will turn off after displaying the last test result.

NOTE

- Any time you wish to exit the memory, keep pressing Main button for 5 seconds or leave it without any action for 3 minutes. The meter will switch off automatically.
- Control solution results are **NOT** included in the day average.
- If using the meter for the first time, "----" displays when you recall the test results or review the average result. It indicates that there is no test result in the memory.

DOWNLOADING RESULTS ONTO A COMPUTER


You can use the meter with a USB cable and the GlucoRx Healthcare Software System to view your test results on your personal computer. To learn more about the GlucoRx Healthcare Software System or to obtain a USB cable separately, please contact GlucoRx customer services or the place of purchase for assistance.

1. Obtaining the required cable and installing the software

To download the GlucoRx Healthcare Software System, please visit GlucoRx website at www.glucorx.co.uk.

2. Connecting to a personal computer

Connect the cable to a USB port on your computer. With the meter switched off, connect the other end of the USB cable to the meter data port. "USB" will appear on the meter display, indicating that the meter is in communication mode.



3. Data transmission

To transmit data, follow the instructions provided with the software. Results will be transmitted with date and time. Remove the cable and the meter will automatically switch off.

Warning:

While the meter is connecting to the PC, it will be unable to perform a blood glucose test.

BATTERY RECHARGE

Your meter comes with one 3V CR2032 lithium battery. If the low battery symbol "⚡" appears on the screen this indicates that the battery is low and it is time to replace the battery.

1. Press the edge of the battery cover and lift it up to remove.
2. Remove the old battery and replace with one 3V CR2032 lithium battery.
3. Close the battery cover. If the battery is inserted correctly, you will hear a "beep" afterwards.

NOTE

- Replacing the battery does not affect the test results stored in the memory.
- As with all small batteries, these batteries should be kept away from small children. If swallowed, promptly seek medical assistance.
- Battery might leak chemicals if unused for a long time. Remove the battery if you are not going to use the device for an extended period (i.e., 3 months or more).
- Properly dispose of the battery according to your local environmental regulations.

DISPLAY MESSAGES AND PROBLEM-SOLVING GUIDE

The following is a summary of display messages. If your meter displays an error message, please follow the actions for the error message as described in the table below. If the problem persists, please contact GlucoRx for help.

MESSAGE	WHAT IT MEANS
Lo	<1.1 mmol/L (20 mg/dL)
Hi	>33.3 mmol/L (600 mg/dL)
Li	1.1 - 33.3 mmol/L (20 - 600 mg/dL)

MESSAGE	WHAT IT MEANS	ACTION
E-b	The symbol appears with this message. This means that the battery does not have enough power for a test.	Replace the battery immediately.
E-U	A used test strip is inserted, or there is an electrical problem.	Repeat test with a new test strip. If the error message re-appears, please contact GlucoRx for help.
E-t	The temperature of the environment, meter, or test strip is outside the system's operating range. You cannot perform a test until the meter and test strip are within the operating range of 10°C to 40°C (50°F to 104°F).	Repeat the test after the meter and test strip are within the operating temperatures.
E-E	Problem with the meter.	Review the instructions and try again with a new test strip. If the problem persists, please contact GlucoRx for help.
E-F	You may have removed the strip after applying blood to the absorbent hole.	Keep the test strip in the meter after it has absorbed the blood until the meter displays the test result.

For questionable results, please see test strip package insert for important information.

SYMBOL INFORMATION

Symbol	Referent	Symbol	Referent
	In vitro diagnostic medical device		Manufacturer
	Consult instructions for use		Serial number
	Temperature limitation		Authorised representative in the European Community
	Use by		
	Batch code		

CARING FOR YOUR METER

Cleaning

1. To clean the meter exterior, wipe it with one disinfecting wipe (Example: Micro-Kill Plus™) to clean exposed surfaces thoroughly and remove any visible dirt, blood or any other body fluid with the wipe. Then dry the device with a soft, dry and clean cloth.

2. Do NOT rinse the meter with water.

SPECIFICATIONS

Model no.: TD-4287

Dimensions: 103.4 (L) x 28.1 (W) x 15.1 (H) mm

Weight: 31.8g

Power source: 3V CR2032 lithium battery

Display: LCD

Memory: 1000 memory sets with date and time

Automatic detection of electrode insertion

Automatic reaction time count-down

Automatic blood detection

Temperature warning

Operating conditions: 10°C to 40°C (50°F to 104°F), below 85% R.H.

Meter storage/transport conditions: -20°C to 60°C (-4°F to 140°F), below 95% R.H.

Strip storage/transport conditions: 2°C to 32°C (35.6°F to 89.6°F), below 85% R.H.

Measurement units: Fixed mmol/L

Measurement range: 1.1–33.3 mmol/L (20–600 mg/dL)

This device has been tested to meet the electrical and safety requirements of: IEC/EN 61010-1, IEC/EN 61010-2-101, EN 61326-1, EN 61326-2-6.

Distributed by GlucoRx Ltd
4 Wintonlea, Monument Way West
Woking, Surrey, GU21 5EN, UK
Website: www.glucorx.co.uk
Customer care number: 01483 755133

TaiDoc Technology Corporation
B1-7F., No. 127, Wugong 2nd Rd., Wugu Dist.,
New Taipei City 248, Taiwan (R.O.C.)
www.taidoc.com

MedNet GmbH
Borkstraße 10, 48163 Münster, Germany



TERMS & CONDITIONS OF WARRANTY

1. We warrant this product to be free of defects in workmanship and materials within the said warranty period on the warranty certificate.
2. During the warranty period, if this product is found to be defective, you may send it, together with the purchase receipt and Warranty Certificate, to our office during normal business hours for warranty service. We will then repair or replace defective parts or exchanging the whole product as we may choose, at no charge to the original owner. After such repair, replacement or exchange, the product will be warranted for up to the remainder of the warranty period.
3. This warranty is valid only if the Warranty Certificate and Warranty Registration Card are duly completed with date of purchase, serial number and dealer's stamp, and if the Warranty Registration Card is sent to our office no later than 6 months from the date of purchase.
4. This warranty is void if this product has been repaired or serviced by unauthorised person. This warranty does not cover defects caused by misuse, abuse, accident, tampering, lack of reasonable care, fire or any other acts beyond human control.
5. Except as stated in the above paragraphs, we disclaim all other warranties, implied or expressed, including the warranties of merchantability or fitness for a particular purpose with respect to the use of this product. We shall not be liable for any direct, consequential or incidental damages arising out of the use or inability to use this product.

PLEASE
AFFIX
STAMP

DEALER'S INFORMATION