

Insulin Pen Needles Universal fit on all Insulin Pens



Glucorx CarePoint pen needles have universal fit on all insulin pen delivery devices including Eli Lilly*, Novo Nordisk AS*, Becton Dickinson*, Owen Mumford* and Sanofi*. They are electropolished with high grade siliconisation for enhanced ease of injection.

These high quality pen needles use Maxflow™ technology (extra thin wall allows maximum flow rate and reduced injection time) and are anti-coring (no fragmentation).

Glucorx CarePoint pen needles have a 3 bevel needle tip.

Glucorx CarePoint



Glucorx CarePoint Ultra



Dear Patient
If you would like to use Glucorx CarePoint Pen Needles, please show this information to your Doctor.

Patient Name

Dear Doctor
Please prescribe me with Glucorx CarePoint Pen Needles. If you require further information, please contact Glucorx on 01483-755133.

☐ Glucorx CarePoint 4mm/32G Ultra Pip Code: 3938966

☐ Glucorx CarePoint 4mm/31G Pip Code: 3759255

☐ Glucorx CarePoint 5mm/31G Pip Code: 3812674

☐ Glucorx CarePoint 6mm/31G Pip Code: 3761582

☐ Glucorx CarePoint 8mm/31G Pip Code: 3761590

☐ Glucorx CarePoint 10mm/29G Pip Code: 3761608

☐ Glucorx CarePoint 12mm/29G Pip Code: 3761681